

OneMinute[®] WB Stripping Enhancing Buffer*

Cat. No. GM6021

1002.25

Number Description

GM6021	<i>OneMinute</i> [®] WB Stripping Enhancing Buffer	25 Tests
	(Buffer A: 25X Supplement 5ml ; Storage: Store at 4 °C)	
	(Buffer B: <i>OneMinute</i> [®] WB Stripping Enhancing Buffer 125ml; Storage: Store at room temperature)	

Please refer to label for the Expiration Date. Product shipped at ambient temperature.

Introduction

Pretreatment of Western blots with *OneMinute*[®] WB Stripping Enhancing Buffer significantly enhances the efficiency of Western Blot stripping reagents (such as *OneMinute*[®] Western Blot Stripping Buffer). Transferred protein on the blot keeps intact after treating with *OneMinute*[®] WB Stripping Enhancing Buffer. The combination of this buffer with *OneMinute*[®] Western Blot Stripping Buffer creates the world-best stripping reagents, especially for high-affinity antibody.

Additional material required

- Western Blot stripping reagents, such as *OneMinute*[®] Western Blot Stripping Buffer (GM6001, GM6002)
- Western blot, previously blocked, probed and detected with a chemiluminescent substrate
- Washing Buffer, such as phosphate-buffered saline (PBS) or Tris-buffered saline (TBS) with 0.05% Tween-20.

Protocol

1. **Wash blot in Washing Buffer to remove the chemiluminescent substrate.** Blot may be stored in PBS or TBS at 4°C until the stripping procedure can be performed. NEVER LET THE BLOT DRY.
2. **Dilute 25X Buffer A into Buffer B to make 1X working solution.** For example, dilute 1ml 25X Buffer A into 24ml Buffer B to give 25ml 1X working solution. 5ml 1X working solution is enough for a standard Western blot (8 x 10 cm blot).
Note: Freshly made 1X working solution is recommended to obtain best results. The 1X working solution is stable for only 2 weeks at room temperature.
3. **Place the blot in 5ml 1X working solution made in step 2, shake for 5min at room temperature. Decant working solution.** Use sufficient volume to ensure the blot is completely wetted by the working solution (5ml per test is required for a standard 8 x 10 cm blot).
4. **Add *OneMinute*[®] Western Blot Stripping Buffer and shake for 30 sec ~ 1min. Then, add 20ml Washing Buffer and decant the mixture.** Use a sufficient volume to ensure the blot is completely wetted by *OneMinute*[®] Western Blot Stripping Buffer (5ml per test is required for a standard 8 x 10 cm blot).
Note: Optimization of stripping condition is essential for best results. The blot can be treated with *OneMinute*[®]

OneMinute[®] WB Stripping Enhancing Buffer^{*}

Cat. No. GM6021

1002.25

Western Blot Stripping Buffer for longer time (~ 5 min) or at higher temperature (37°C).

5. **Wash the blot in 10~20ml Washing Buffer for 3 times.** Briefly shake in hand for 10~20 seconds each time.
6. **The blot is ready for subsequent re-probing.** Perform the re-probing with primary antibody in proper buffer containing 1~5% milk. Re-blocking the blot is usually not necessary but might be required in some cases.

Notes:

- Be sure to keep the Western blot wet all the time. Never let the blot dry.
- Freshly made 1X working solution is recommended to obtain best results. The 1X working solution is stable for only 2 weeks at room temperature.
- Stripping and re-probing fluorescent Western blots is not recommended because unexpected background is generated.

This product is for research use only and not for human or diagnosis use. The warranty limits our liability only to replacement of the product. There is no obligation to replace this product as the result of accident, disaster, misuse, fault, negligence, improper storage or handling of the products. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by GM Biosciences Inc. We shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use the product.

* US Patent Pending